

## IMPORTANT PRODUCT INFORMATION

### Resistance against the cleaner „Bleach“

DEAR CUSTOMER,

The products shown below are tested successfully in our test laboratory for resistance against solutions out of household bleach (with max 50,000 ppm of free chlorine) in a dilution of 1:5. A Solution with this high concentration of 10,000 ppm of free chlorine, exceed valid requirements of 6,000 ppm of CDC (Centers for Disease Control and Prevention)-Guidelines for Biosafety in Microbiological and Biomedical Laboratories 5<sup>th</sup> Edition.

„Bleach“ is a chemical that whitens clothing. It is also used as household chemical to remove stains and as disinfectant, primarily in the bathroom and kitchen.

INDUPROOF 3, BLACK



TKG-086-MB-IP68-BLACK

INDUPROOF MED, GREY



TKG-105-MED-IP68-GREY

INDUPROOF MED, BLACK



TKG-105-MED-IP68-BLACK

SMART CLINICAL BOARD



TKG-106-IP68-WHITE

SMART CLINICAL BOARD - COMPACT



TKG-086-IP68-WHITE

SMART CLINICAL BOARD - TOUCH



TKG-107-TOUCH-IP68-WHITE

INDUPROOF PRO - MAG



TKG-107-TOUCH-IP68-WHITE-MAG

INDUPROOF PRO - MAG + BACKL  
- BACKL




TKG-107-TOUCH-IP68-WHITE-(MAG)-BACKL

INDUMOUSE, BLACK / WHITE / GREY




TKH-MOUSE-IP68-SCROLL-OPT

CLEANKEYS CK3



CK3PGW

CAPACITIVE KEYBOARD



TKR-103-TOUCH-ADH(C)

COMPACT CAPACITIVE KEYBOARD



TKR-086-TOUCH-ADH(C)